



Kits & Accessories

Clinical Electrophoresis

Details Info

Features:

Unlike dry cellulose acetate strips that are restricted to low volume micro tests, CellasGEL's greater thickness allows tests to be performed on semi-micro and macro scales using a wide range of specialist applicators. Consequently greater sample-volumes may be loaded as larger, but finer bands over a wider front. This reduces sample saturation and aids densitometric band quantitation, thereby improving resolution.

Ordering Information:

ACAPPS22	Applicator Semimicro 2+2 for CellasGEL 2.5x14 cm and Kit code ACKITCU
ACAPPS4SP	Applicator Semimicro 1 X 4 specimen semi-micro applicator
ACAPPS4	Applicator Semimicro 4 for 1.2 _l/9mm – First type
ACAPPS6	Applicator Semimicro 1 X 6 specimen semi-micro applicator
ACAPPM8	Applicator Micro 8 for CellasGEL
ACDENS	TurboScan Software Densitometer (computer and scanner not included)
ACSCAN	Flatbed Scanner for TurboScan software.
ACBDG8.5S	ACBDG8.5S Bridge 8.5 cm (standard size for CellasGEL kits) for CellasGEL strips 5.7x14 or 2.5x14cm
AC-CAB	Electrophoresis cable (Black & Red)

Specification:

ACKITCU	CellasKIT IFE SERUM + CONCENTRATED URINE (5+5 tests for 5 patients) on CellasGEL
ACKITI2432	ACKITI2432 CellasKIT IMMUNOFIXATION (24 semimicro or 32 micro tests) on CellasGEL
ACKITSP100200	CellasKIT SERUM PROTEINS (100 semimicro or 200 micro tests) on CellasGEL
ACKITSP150200	CellasKIT HIGH RESOLUTION OF SERUM PROTEINS (150 semimicro or 200 micro tests) on CellasGEL
ACKITHG100	CellasKIT HEMOGLOBINS (100 semimicro tests) on CellasGEL
ACKITLP100	CellasKIT LIPOPROTEINS (100 semimicro tests) on CellasGEL
ACKITCU-ABS	Includes 1 x ACKITCU, 1 x ACAPPS22 & 3 x ACBDG8.5S
ACKITI2432-ABS	Includes 1 x ACKITI2432, 2 X ACAPPS6 & 6 X ACBDG8.5S
ACKITI2432-ABM	Includes 1 x ACKITI2432, 2 X ACAPPM8 & 6 X ACBDG8.5S
ACKITSP-ABS	Includes 1 x ACKITSP100200, 1 X ACAPPS4 & 3 X ACBDG8.5S
ACKITSP-ABM	Includes 1 x ACKITSP100200, 1 X ACAPPM8, 3 X ACBDG8.5S
ACKITSPHR-ABS	Includes 1 x ACKITSP150200, 1 X ACAPPS4, 3 X ACBDG8.5S
ACKITSPHR-ABM	Includes 1 x ACKITSP150200, 1 X ACAPPM8, 3 X ACBDG8.5S
ACKITHG100-ABS	Includes 1 x ACKITHG100, 1 X ACAPPS4, 3 X ACBDG8.5S
ACKITLP100-ABS	Includes 1 x ACKITLP100, 1 X ACAPPS4, 3 X ACBDG8.5S



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: