



UAZ20

| Urine Analyzer

☰ Details Info

Specification:

| | |
|--|--|
| Test item | WBC, urobilinogen, bilirubin, specific gravity, glucose, PH, protein, occult blood, nitrite, ketone bodies, vitamin C, Albumin |
| Strip type | 8, 10, 11, 12 optional, strips open |
| Apply strip | URS-10A, URS-11A, F11-II, F12 |
| Test speed | 120 / hour, 60 / hour available |
| Data Storage | 1000 results |
| Display Screen | 240 x 64 LCD screen |
| Automatic residue suction system to prevent cross-contamination between samples. | |

Automatic fault detection and alarm system to ensure accurate measurement results.
Automatic correction of nonspecific test strip, urine pH, specific gravity and color have
the impact on the test results

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved