



AC-Chiller

Comet Assay Tanks

Details Info

Features:

The AC-CHILLER is available ready assembled with the thermostat mounted on the refrigerator and supplied with insulated tubing and clips to form a system ready to use. Supported by an industry leading 3 year warranty emphasising our confidence in the reliability of our product.

The system comprises of the following:

- 5 litre tank with drain tap
- relay control for refrigeration on/off
- cooling power 140W @ 0°C
- temperature range -20 to 100°C

Model	AC-CHILLER	Chiller System, -20 - 100°C, inc tubing & clips - 240V
	AC-CHILLER\$	Chiller System, -20 - 100°C, inc tubing & clips - 115V
	F250	Julabo Chiller

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com