



DCU-3190

| Dental Unit

Details Info

Features:

- ◆Humanized design concept adopts the theory and method of 3D simulation and ergonomic design.Combined with the person and machine engineering theory and three-dimension dynamic simulation design.
- ◆ Advanced digital technology which is good for high efficiency and reach the accurate control for all the instruments.
- ◆ Health and environmental protection concept: Automatic pipe sterilization system and anti-reaction.
- ◆ Combine with high quality spare part and advanced produce process, the dental unit adopt innovation idea and modularization design.
- ◆ Strong LCD function for dialogue between person and machine.
- ◆ Practical water shortage warning system.
- ◆ Instrument tray can be enlarged by needs.

Standard accessories			
	Electriccontrol	Air control	Separated
Fibre optic handpiece	1pcs	1pcs	1pcs
Low speed electric motor	1set		
Newtron Scaling Machine	1set		

Import LED Curing Light	1set		
FARO Operation Light	1set	1set	
LED Sensor Operation Light			1set
3 way syringe	2set	2set	2set
Clean water system	1set	1set	1set
Automatic warm water& rinsing system	1set	1set	1set
Motor driven dental chair	1set	1set	1set
Clinic stool	1set	1set	1set

Packing size	Down
Size (mm)	1840X1230X112
) kg (N/W	245
) kg (G/W	365

Notes:

- ◆ 16 illumination leves for optical handpiece.
- ◆ Water acquisition by sensor or manual.
- ◆Chair safety system for dental chair's movement and assistant stand which is convenient for operation.
- ◆ Hand holder is convenient for 4 hand's operation.
- ◆Suction delay function.
- ◆Multi-functional foot control
- ◆ Operation light control.
- ◆Free revolution of telescope.
- ◆ Movement of dental chair can be controlled.
- ◆ Multi-light source realizes the shadow-free operation.
- ◆ Tri-axial rotation design ensures the complete illumination.
- ◆ Sensor stepless adjustment makes long lifespa.
- ◆ The angle of handle can be adjusted freely, and also detachable and easily disinfected.
- ◆ Pipe disinfection system avoid cross infection.



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: