



## DCU-3195

| Dental Unit

☰ Details Info

DCU-3195 - Implant (standard)

Features:

DCU-3195 Dental unit is use for implant, operation, which meets implant and operation's needs. Workbench, Electric surgery dental bed, shadow less operating light, be sterilized suction system, portable assistant workbench are included to meet the dentist's needs.

- ◆ Soft star smooth movement random vibration.
- ◆ Pear-shaped backrest fits for surgery in different way. Workbench for assistant
- ◆ Stainless steel material, which easy for disinfection.
- ◆ Intermittent work, Interlock Switch, foot control.
- ◆ Suction, saliva system and handling funnel cup.
- ◆ Suction ,saliva and handling spittoon.
- ◆ No discharge system, but two 2.5L waste bottle installed on device.
- ◆ Moveable and locked instrument cart. Working bench for dentist
- ◆ Stainless steel material, which easy for disinfection.
- ◆ Electric lifting instrument ( height cope between 730mm~850mm ).

- ◆ 2 high speed hand pieces, 1 low speed handpiece, 1 three way syringe.
- ◆ Moveable and locked instrument cart.
- ◆ Easy for any brand of Implant unit's installation. Shadow less operating light
- ◆ 5 sections for LED lamp, intensity is 50000lx.
- ◆ Different shadow less lamp fixed height meet consumers 'clinic installation needs shadow less operating light from ceiling or dental unit mounted

Standard accessories	
3 way syringe	1pc
Lamp with column	1set
Working bench for dentist	1set
Workbench for assistant	1set
Suction system	1set
Working bench for dentist	1set
Clinic stool	1set



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved