



DCU-3200

| Dental Unit

☰ Details Info

Features:

- ◆Wireless foot control, carrying out linearly adjustment of the moveable instrument
- ◆Dentist, assistant position integrate the 7" TFT color touchable screen
- ◆Memory procedures' supply for 3 dentists' individual set
- ◆Soft-start for Dental chair
- ◆Oral disease adopts remote control
- ◆Professional remote teaching for oral medicine
- ◆Advanced self-warning and inspection by remote control system

Standard accessories	
Fibre optic handpiece	1pc
Electric motor	1set

Import air motor	1set
Curing Light	1set
FARO Operation Light	1set
X film viewer	1set
Three way syringe	2pcs
Pipe sterilization system	1set
Automatic warm water& rinsing system	1set
Digital oral camera	1set
Motor driven dental chair	1set
Clinic stool	1set

Packing size		
Item	(mm)	gk
Dental unit	1600×780×1200	180
Dental chair	1700×1200×1160	320

Notes:**Automation:**

- ◆Electric lifting system for headrest
- ◆Automatic Adjustment for dental unit
- ◆Electric Split Spittoon

Intelligentization:

- ◆Untouchable controlling table patent
- ◆Multiple chair mounted Unit (light mounted chair, cup mounted spittoon, spittoon mounted with rising position)
- ◆Dual sensor rinsing water function

Avoid cross-infection:

- ◆Pipeline anti-suction system
- ◆Pipeline disinfection system
- ◆Pipe adopts anti-retraction system
- ◆Pipe sterilization system
- ◆Non-contact manipulation
- ◆Suction delay functions
- ◆Untouchable controlling system
- ◆Time-expanding for saliva & suction system



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

