



Dot and Slot Blotters

| Blotting

Details Info

Four different sample number size and style of Hybridisation Manifold are offered:- two types of dot blotter and two types of slot blotter.

Typical applications include clone screening with DNA / RNA probes in Southern/northern blots and immunological screening with antibodies in western blots.

The units incorporate precision lapped surfaces to ensure uniform blotting membrane contact and a leak proof gasket. These prevent lateral transfer of samples – smudging – by ensuring that a complete vacuum is formed. Six thumbscrews ensure even and tight sealing for fast sample suck down.

Dot blotters are available in 48 and 96 well versions and slot blotters in 24 and 48 well.

AC-D48	Dot Botter, 48 Sample
AC-D96	Dot Blotter, 96 Sample
AC-S24	Slot Blotter, 24 Sample
AC-S48	Slot Blotter, 48 Sample

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com