



MRM35- Rotary Paraffin

| Microtome Rotary

Details Info

Features

- 1. Ergonomically one-step Injection Molded housing results in appropriate structure, attractive appearance and convenient operation.
- 2. Large magnetic waste tray enables easy removal for cleaning, and was designed black to blend into one harmonious whole together with the main machine.
- 3. Hand-wheel support multi-locking evenly around the wheel
- 4. The blade clamp can be moved horizontally to avoid danger from unintended touch on blade.
- 5. Guide rails made in USA are so highly precise that enables to slide with even no gap, thus results in precision, longer equipment life, free of noise.
- 6. Fast rewind function.
- 7. Indicator displays clearly scales

Specifications:

Section thicknessrange	1~60μm(3~5μm thickness leads to best section results)
Section thickness setting	0.5~2μm,in 0.5μm increment 2~10μm,in 1μm increments 10~20μm,in 2μm increments 20~60μm,in 5μm increments
Specimen vertical stroke	35mm

Specimen horizontal stroke	23mm
Max specimen size	60mm×50mm×40mm
Specimen positioning	horizontally 8°, vertically 8
Dimensions (W×L×H)	55mm×46mm×428mm
Weight	About 25kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N  
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

[sales@acculabusa.com](mailto:sales@acculabusa.com) [middleeast@acculabusa.com](mailto:middleeast@acculabusa.com)



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved