



MRM40- Manual Rotary

| Microtome Rotary

Details Info

Features

Spacious magnetized section wastr tray ,Ergo-handwheel handle

Convenient storage area on top of instruments ,Disposable blade holder with intergraded satety finger protection guard

Two handwheel locks :handwheel can be locked at any position,Distortion – resistant base plate ensures optimum overall stability,Both low&high profile disposable blade are usable

Two types of specimen clamp ;standard specimen clamp and “C” type specimen clamp

Retraction function :220μm,could be turned on \off

Extremely smooth running force compensation system

Specification:

Section thicknes setting	0 μm ~ 60 μm adjustable 0 μm ~ 2 μm ,in 0.5 μm increments 2 μm ~ 10 μm ,in 1 μm increments 10 μm ~ 20 μm ,in 2 μm increments 20 μm ~ 60 μm ,in 5 μm increments
Trimming thickness	10 ~ 50 μm adjustable
Retraction range	220 μm , could be turned off
Maximum specimen size	50 x60 x40 mm or standard cassette

Specimen vertical stroke	70 mm
Specimen horizontal stroke	30 mm
Specimen holder adjusted system	horizontal orientation : $\pm 8^\circ$ vertical orientation : $\pm 8^\circ$
Blades holder base reposition (left – right)	50 mm
Dimensions	width :460mm (320mm withouthandwheel) Depth :550mm Height ;300mm
Weight	about 42 KG
1unit/1 plywood case size	56*48*46cm
N.W:42kgs G.W: 50kgs	



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: