



Mult Horizontal Screen

| Horizontal Electrophoresis

☰ Details Info

Product Description

Designed for rapid screening of very large numbers of Clonal or PCR samples, the Multi Horizontal Screen horizontal gel unit has a 672 maximum sample capacity per gel. This allows loading and analysis of seven 96-well format microtitre plates. The large gel run length of 32cm also allows resolution of samples over a long distance for separation of complex sample bands, such as in restriction fragment analysis.

The unit is available with a full length tray or with other tray length options of 16 or 24cm so that the user's exact requirements can be matched. In addition to options for single length gel trays, the MultiHorizontal Screen is available with all three gel tray lengths to provide the maximum in flexibility, versatility and value. Buffer recirculation ports are included as standard to allow enhanced resolution over extended runs while loading guides enhance well visibility for easy sample loading.

Specifications

Gel dimensions (w x l)	26 x 32cm, 26 x 24cm, 26 x 16cm
Unit dimensions (w x l x h)	28 x 50 x 9cm
Max. sample capacity	26x32cm tray - 672 samples, 26x24cm tray - 504 samples, 26x16cm tray - 336 samples
Buffer volume	1400ml

Combs available:

No. of samples	28MC, 56MC
Thicknesses	0.75, 1, 1.5, 2mm

Ordering information :

Mult Horizontal Screen	
ACSCREEN32	MULTI HORIZONTAL Screen, 26 x 32cm UV Tray, 6 x 28 sample combs, loading guides & MS26-FC
ACSCREEN24	MULTI HORIZONTAL Screen, 26 x 24cm UV Tray, 6 x 28 sample combs, loading guides & MS26-FC
ACSCREEN16	MULTI HORIZONTAL Screen, 26 x 16cm UV Tray, 6 x 28 sample combs, loading guides, & MS26-FC
ACSCREENTRIO	MULTI HORIZONTAL Screen, 26x16cm, 26x24cm, 26x32cm UV Trays, 6 x 28 sample combs, loading guides & MS26-FC
AC26-UV32	26 x 32cm UV Tray
AC26-UV24	26 x 24cm UV Tray
AC26-UV16	26 x 16cm UV Tray
AC26-LG	Loading Guides
AC26-WP	Viewing Platform
AC26-PE	Positive Electrode
AC26-NE	Negative Electrode
AACCRNCP	Cool-pack and platform for MSSCREEN
AACCRNBACSB	Buffer Saver Blocks for MSSCREEN, pk/2 saves 625ml of buffer
AC26-FC	Flexicaster for gels up to 32cm. Casts 7, 10, 15, 16, 20, 24 and 32cm long gels
AC26-UACV	26cm UV Gel Scoop
AC26-28MC-0.75	Comb 28 sample MC, 0.75mm thick
AC26-56MC-0.75	Comb 56 sample MC, 0.75mm thick
AC26-28MC-1	Comb 28 sample MC, 1mm thick
AC26-56MC-1	Comb 56 sample MC, 1mm thick

AC26-28MC-1.5	Comb 28 sample MC, 1.5mm thick
AC26-56MC-1.5	Comb 56 sample MC, 1.5mm thick
AC26-28MC-2	Comb 28 sample MC, 2mm thick
AC26-56MC-2	Comb 56 sample MC, 2mm thick
AC26LID	MS26 lid (No cables)
AC26TANK	MS26 Tank (including electrodes)
AC26-ST	Gel tray sealing tape. 65m X 25.4mm
AC-CAB	Electrophoresis cable (Black & Red)



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved