



## PDU-3020ST

| Dental Unit

Details Info

### Dental unit(Standard-Model)

Features;

PDU-3020ST Dental unit (Standard-Model) is a fully self-contained electric dental system.

This type of dental unit adopts high quality alloy aluminum materials and including oil less compressor, tank, water-air-electric control system. Quick and easy to set up and operate. Ready for travel It is a section of economical practical dental unit.

Volume: 510x360x720mm

Weight: 15-25kg

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com