



PFT 29

| Paraffin Trimmer

☰ Details Info

Features:

- Innovative design help user to remove excess paraffin from cassettes easily , quickly and safely
- Remove the excess paraffin by scraping the block on its heated surface .
- Its temperature could be set from ambient to 90°C.
- It is easy and safe to be opearted by a user friendly interface .

Specification:

Supply voltage	220V±10%, 50hz /110V±10%, 60hz
Table top (mm)	150×220
Temperature control accuracy	±1°C
Outside dimensions (W*D*H)	220×160×175 mm
Weight	about 2.5 KGS



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

