



PSC6 Paraffin stone cabinet ( 6 drawers with lock )

| Slide Storage Cabinet

Details Info

Specification information :

The base use Cold-rolled steel plate SPCC1.2 from Baosteel; the cabinet use Cold-rolled steel plate 0.8; the drawer use Cold-rolled steel plate 0.8.  
The biomedical paraffin block cabinet has three independent cabinet, each cabinet has six drawers, totally 18 drawers. Paraffin wax block can be put freely with its size. It can store totally around 11,000 blocks.

|                    |   |
|--------------------|---|
| Size               | H1295 x W450 x D478 mm  |
| Weight             | 80 kgs  |
| Surface treatment: | Defatted oil, Table regulator, Zinc phosphating, passivation, Powder spraying.. |
| Slot:              | The slot is made by special mold  |
| Mark groove        | Integration forming   |
|                    |   |

Specialized slide drawer, ABS slide plate



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N  
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved