



PSP PATHOLOGY SAMPLING SECTION

| Pathological sampling bench

Details Info

Advantage

- 1. Odor-remove function.
- 2. Single water sink, two-person operations.
- 3. Exhaust air system: High-power suction device, low noise.
- 4. The base cabinet is inward-concave type, convenient for sitting operation.
- 5. Automatic flushing device, slightly right-side inclined worktable, easy for waste discharge and cleaning.
- 6. Optional water heater is available, realizing hot &cold water adjustment.
- 7. Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes.

Specification:

Model	PSP-150	PSP-180
External Size (W*D*H)	1500*750*1950mm	1800*750*1950mm

Working Area Size (W*D*H)		1395*430*900mm	1695*430*900mm
Work Surface Height		900mm	
Consumption		300W	
Power Supply		110V/220V±10%, 50/60HZ	
Exhaust Airflow Volume		500-600m3/h	
Material	Main Body	304 stainless steel, thickness 1.2mm	
	Work Table	304 stainless steel, thickness 1.5mm	
Standard Accessories		Spotlight UV lamp Fluorescent lamp Exhaust system Auto flushing device Stainless steel hanger Adjustable hot and cold water tap 180mm deep stainless steel water sink, funnel with filter.	
Optional Accessories		Pulverizer, Electrical water heater, Automatic induction hand dryer, Eye washer, Camera.	



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
 Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: