



SDT-Sries Smart disintegration tester

| Pharmaceutical industry

Details Info

Specification:

Intelligent Control Controls alarm	
A clock timer function	
Take the next heating water bath	
Overall neat appearance	
Instruments owned lighting device	
Temperature-controlled tanks	37 °C ± 1 °C

Gondola lift frequency	30-32 beats / min
Pharmacopoeia	2010 edition of the bezel meets the provisions of
Power	220v / 50Hz / 500w



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved