



SST18

| STAINER

Details Info

Specification			
1	Quantity of tank	18 PCS	First tank is cleaning tank (flowing water) The 18 th tank is entering tank (heating). The 17 th tank is ending tank.
2	Capacity	1000ml	can be set 47 slides
3	Power supply	220V±10% 50Hz。	
4	Power	Max. 400W	
5	Set temperature	0-99℃	
6	Run mode	Run the baskets continuous	Can stain large batch slides continuous
7	Retention time in each tank	1 second to 99 Min. and 59 Sec.	

8	Store Program	6 sets	
9	Dimension	140×50×55cm	
Function			
1	Run baskets continuous	Can stain large batch slides continuous	
2	Colorized touch screen	Display all running parameters, easy to watch and convenient to operation.	
3	With ion filter		
4	Store programs	Can store 6 sets programs	
5	Movement parts	Movements parts choose bearing linear guide, high precision, small resistance, steady running, strong abrasability.	
6	Water washing	Flowing water to wash the basket, fill and drain water automatically	
7	With heating, for drying slides function		



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: