



TBP20

| Tissue Bath & Dryer

Details Info

Feature :

- Adopted the newest high technology and nice screen control, new type of heating material to ensure Small volume and light weight.
- Automatic memory of setting temperature and running .
- High precision in temperature controlling. heating fast ,durable and energy-serving .
- If the power is switched off suddenly, the machine can go on work and correct mistake automatically when the power switch on again.
- The LCD display the current temperature and presetting temperature to make the communication between user and machine easily .

Technical specification:

Temperature prefabricate	0~100℃ set arbitrarily
Rating Power	300W
Power voltage	AC220V±10 % 50HZ
Dimension	425×415×215mm
Thermostatic control with an accuracy of ±0.5℃	

Weight	10kg
--------	------



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: