



TFT- Series Tablet friability tester

| Pharmaceutical industry

Details Info

Tablet friability tester

Tablet friability tester is consist of the control system, transmission, wheel parts and other components, precision control system by using singlechip microcomputer composed of parts centralized control instruments reasonable structure, high degree of automation, control progress, high sensitivity, easy to operate, reliable work.

type	TFT-300	TFT-100
Cylinder internal dimensions	39mm diameter 286mm depth	39mm diameter 286mm depth
Fall height	156mm	156mm
Speed Accuracy	±0.1 rpm	±0.1 rpm

Speed Time	Adjustable 0-300 rev 9999 min	Fixed type
Working hours	100 rpm 4 min 25 r/min	100 rpm 4 min 25r/min
Dimensions	360*300*340mm	360*300*340mm



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office  
: warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

[sales@acculabusa.com](mailto:sales@acculabusa.com) [middleeast@acculabusa.com](mailto:middleeast@acculabusa.com)



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved