



THT-100 Tablet hardness

| Pharmaceutical industry

Details Info

Features:

300 Tablet hardness tester is the special test equipment to test the hardness of a tablet drug quality , this instrument is suitable tablet hardness testing of drugs, are pharmaceutical, medical teaching and research, laboratory research and drug sector of essential equipment.

The instrument uses a high-speed single-chip microcomputer control system, there are powerful data processing capability, automatic test maximum tablets, minimum , average, standard deviation and coefficient of variation reports. The instrument has a user-friendly blue LCD large screen graphic LCD display, operating functions using the menu mode, easy to operate, truly interactive content and clear, high precision, test repeatability, intuitive readings clear, the unit of measurement : kg-N can be set in the software operation.

Specification:

Tablets pressure measuring range	1 to 30 kg
Determination of the maximum diameter of the tablet	30mm

Hardness Pressure accuracy	± 0.5% of maximum range
0.01kg	0.01kg
Input power	30W



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: