



WBG10X10 Western Blotting

| Blotting

Details Info

Features:

WBG10X10 Trans-Blot is compatible with the VGE10X10D vertical electrophoresis tanks meeting the general need of protein and Nucleic tests. With special design and injection moulding technology, two gels can be run simultaneously under identical buffer conditions with high-resolution separation and simple operation.

The WBG10X10 Trans-Blot adopts a safety design concept with automatic power cut-off when the lid is removed.

Specification:

Gel transfer size(W×L)mm	100×100.
Buffer Volume	750ml

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com