



WBG17X16 Western Blotting

| Blotting

Details Info

Feature:

WBG17X16 Trans-Blot is generally intend for the transfer print when the run is finished.
The Protein or Nucleic Acid is moved to Acetate Cellulose by using the Trans-blot system.

Specification:

Gel transfer size	(W×L)mm:170×160.
Buffer Volume	4000ml

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:
sales@acculabusa.com middleeast@acculabusa.com

Phone Support: