



## ACZ-152 Single-channel

| Coagulation Analyzer

☰ Details Info

ACZ-152 Single-channel

### Features:

Testing channels:1 channel  
Sample pre-warming position: 5 positions  
Reagent pre-warming position: 2 positions  
Advanced scattered light principle and percentage analysis ensure accurate results.  
Low reagent consumption less than 20uL, open reagent.  
Special testing cup location system, adding sample bracket  
Advanced optical light matching system ensure precision result.  
Optional PT Derived method for testing fibrinogen, decrease the dosage of reagent  
Top quality thermostat, precision 37 °C, ± 0.3  
Dispense precision: Cv<2%  
Memory for 10.000 testing results  
Internal thermal-sensitive printer  
Electronically-linked pipette  
Perform all routine tests such as PT, APTT, TT Fibrinogen and clotting factors  
RS232 interface, PC connecting