



## DNA Sequencing Units

DNA Sequencing Units

Details Info

### Features:

- ACQ Series Nucleic acid Electrophoresis Sequencing Cells are designed for AFLP analysis, difference displays and SSCP research and also fit for heteroduplex, DNA foot printing, RNase protection test, S1 Nucleic acid anzyme figures and primer extension experiments.
- As an experienced manufacturer, acculab has established self-owned brand and models.
- ACQ Series Nucleic Acid sequencing products, including 3 models: ACQ20 , ACQ33, ACQ35.
1. Simple operation ,only four fastening knobs needed
  2. 50% buffer savings compared with traditional tanks
  3. Unfolded electrodes make repair and parts replacement easy
  4. Special heat dissipation processed aluminum plate firmly clings to the gel plates , which effectively guarantees
  5. the gel temperature is even , ensuring the electrophoresis trip is straighter
  6. Auto -switch -off when the lid is removed

| Model                  | ACQ20        | ACQ33       | ACQ35       |
|------------------------|--------------|-------------|-------------|
| Glass plate area (WxL) | 200X450 (mm) | (mm)350X450 | (mm)350X350 |

|                 |             |              |              |
|-----------------|-------------|--------------|--------------|
| Gel area (WXL)  | 180X420 mm  | (mm)330X420  | (mm)330X320  |
| Gel thickness   | (mm)0.4     | (mm)0.4      | (mm)0.4      |
| Amount of gel   | 1(piece)    | 1(piece)     | 1(piece)     |
| No.of samples   | 32.39       | 50,62,74,100 | 50,62,74,100 |
| Buffer capacity | 500(ml)     | 750(ml)      | 750(ml)      |
| Volume (LxWxH)  | 310X230X500 | 455X230X500  | 460x230x400  |
| Weight          | 11.0(kg)    | 16.5(kg)     | 13.2 (kg)    |



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

[sales@acculabusa.com](mailto:sales@acculabusa.com) [middleeast@acculabusa.com](mailto:middleeast@acculabusa.com)



Phone Support: