



ESR-100

| ESR Analyzer

Details Info

Features:

- ESR (westergren and wintrobe Value) and HCT
- ESR data reading in 60 minutes and 30minutes
- 100 channels with random independent test
- Colorful LCD with touch screen
- Barcode scanner and LIS/HIS connection
- Win CE system with unlimited data storage

Specifications:

Measuring Principle	By digital infrared Kinetic photometry
Test Channels	100
Measuring Items	ESR (wintergreen and winrobe Value) HCT
ESR Measuring Range	(0-200)mm/h
ESR Measuring Accuracy	±1mm
HCT Measuring Range	0.2~1
HCT Measuring Accuracy	±0.01

CV	≤3%
Sample Volume	1.5ml
ESR Measuring Time	30 minute or 60 minute selectable
HCT Measuring Time	<1 minute
Date Inquiry	By items (index, data and so on)
Loading pattern	Random
Active Pattern	Automatically detect and activate
Input	Colorful LCD with touch screen
Data Storage	Unlimited with Flash Disk
Screen	Colorful LCD
Display	ESR real time result and curve on LCD
Printer	Internal thermal printer
Printing Content	Result& curve
Data port	RS-232 / USB / LAN port supporting HIS/LIS
Data uploading	Auto and Manual
OS	Android



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved