



Electroblotting transfer Buffers

Reagents & Chemicals

Details Info

Towbin Buffer Tris-Glycine base is available in powder sachets to make 1 litre of 1x working solution for wet and semi-dry transfer techniques. Simply mix 800ml of 1x working solution with 200ml of methanol for complete 1x Towbin Transfer Buffer is available as a 10x concentrate and may be mixed with 900ml of distilled water to make a 1x working solution.

AC-TTGP	Powdered Towbin Buffer, Tris-Glycine base - 10 pouches/pk(10liters/pk)
AC-WTB	10x Western Transfer Buffer (100ml)

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE
Email Us At:
sales@acculabusa.com middleeast@acculabusa.com
Phone Support: