



PS Series Universal

| Power supply

Details Info

Features:

Molding shell ,touch keys ,dual core microprocessor intelligent control:
Indicate the preset value and the actual output value at the same time
It can store 10 electrophoresis method
With automatic memeory function
With standard ,voltage ,constant current ,constant power,
Misoperation ,fault intelligent prompt functions
Automatic detection of no-load ,short cicuit, rapid resistance change , ground leak and system overheating

Specification:

Model	PS2204	PS344	PS3154	PS654	PS614
Voltage range	5-200V	3-300V	5-300V	10-600V	10-600V
Current range	1~2000 mA	1~400 mA	1~1500 mA	1~500 mA	1~1000 mA
Power range	1~200W	1~20W	1~450W	1~300W	1~300W

Type of control	constant-V , Constant I and Constant P	
Display type	LCD screen with backlight	
Output jack	Four sets	
Storage function	Store 10 method	
Timing function	1 min -99h59min	
Power recovery function	Optional	
Automatic memory function	yes	
Automatic shutdown function	yes	
Intelligent prompt functions	yes	
The non –proliferation field function	Optional	
Safety performance	yes	
Volum(LXWXH)	280X237X118	263X275X125(mm)
Weight	3.2kg	4.4kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: