



UV Sterilisation Cabinet

UV Sterilisation Cabinet

Details Info

These UV Sterilization Cabinets provide a convenient area for setting up PCR reactions in a nucleic acid free environment thus limiting PCR reaction contamination. Acting effectively as a low cost alternative to a clean room, the powerful UV lights on the cabinets denature nucleic acids in 5 to 30 minutes making them unsuitable for amplification. The cabinet incorporates safety features to prevent user exposure to UV light. The UV lights themselves are timer controlled and there are safety switches on the cabinets doors which turn off the UV lights when opened. The cabinets also act as an efficient shield of beta radioactive emissions and can therefore be safely used with isotopes such as ³²P. White light provides excellent visibility when working within the cabinets. For use with trays of size 68 x 54cm, CSR-TW4/TY4/TL4. With the benefits comparable to those offered by the AC-UVCAB, the new Mini UV Cabinet is the perfect solution for laboratories constrained by space and budget. A reduced height and footprint allows the Mini UV Cabinet to fit snugly and conveniently into those small areas within the laboratory that would otherwise remain unused.

Ordering information:

AC-UVCAB	UV sterilisation cabinet with timer, four UV lights and white light, no Tray - 230V
AC-UVCAB\$	UV sterilisation cabinet with timer, four UV lights and white light, no Tray - 110V
AC-UVCABTY4	UV sterilisation cabinet with timer, four UV lights, white light, CSR-TY4 Tray - 230V
AC-UVCABTY4\$	UV sterilisation cabinet with timer, four UV lights, white light, CSR-TY4 Tray - 110V
AC-UVCABMINI	Mini UV sterilisation cabinet with timer, four UV lights and white light, no Tray - 230V
AC-UVCABMINI\$	Mini UV sterilisation cabinet with timer, four UV lights and white light, no Tray - 110V

AC-UVCABMTY4	Mini UV sterilisation cabinet with timer, four UV lights, white light, CSR-TY4 Tray - 230V
AC-UVCABMTY4\$	Mini UV sterilisation cabinet with timer, four UV lights, white light, CSR-TY4 Tray - 110V

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support: